

## TECHNICAL DATA SHEET

### InterMask™ : KN95 Face Mask, 12889-001, Protective Mask

#### Intended Use

InterMask KN95 Face Mask, Disposable Non-Woven are intended to be wear it in the daily operation of ordinary environment, covering the mouth, nose and jaw of the user, in order to prevent dust, haze, direct penetration of pathogenic microorganisms, body fluids, and particulate matter provides a physical barrier. The highly efficient filtration media can filter or trap all the common bacteria and particle surrounding us.


12889-001 KN95 Face Mask, Protective Mask			
Characteristic		Results	
Breathing Resistance		Meets FFP2	
Exhalation Valve		No	
Formaldehyde Content(mg/kg)		Not Detected	
Particle Filtration Efficiency (PFE)		>99%	
pH Value		6.5	
Number of Layers		5	
Characteristic	Specification		
Materials	<b>Mask Body Composition</b>	Outer Layer 1	Non-Woven Fabric
		Outer Layer 2	Non-Woven fabric
		Inner Layer 1	Hot Air Cotton
		Inner Layer 2	Melt Blown Woven Fabric
		Inner Layer 3	Non-Woven Fabric
	<b>Nose-Barrett</b>	Adjustable	
	<b>Earloop</b>	Elastic non-woven fabric	
Description	Particle filtering half mask covering the nose, mouth & chin, designed to protect against both solid & liquid aerosols.		
Earloop Style	Elastic non-woven fabric		
Dimensions	<b>Body Size</b>	Length: 234 mm (+/- 1cm)	
		Width: 136.7 mm (+/- 1cm)	

## TECHNICAL DATA SHEET

### InterMask™ : KN95 Face Mask, 12889-001, Protective Mask

Length of Earloop	Elongation between 220% and 260%
Pleat Depth	

#### Visuals



#### Additional Information

Storage Condition	Factory Certifications	Country of Origin
Store under cool, clean, and dry conditions. Avoid excessive heat (over 40°C or 104°F)	GB2626-2006 / FFP2	China

Features	
Outer Layer 1 & 2	Polypropylene Non-woven fabrics, Water-repellent, to prevent fluids from entering the mask
Inner Layer 1	Electrostatic Adsorption Filtering Cotton.
Inner Layer 2	Soft and comfortable PP non-woven
Inner Layer 3	Skin-friendly PP Non-Woven Fabric

## TECHNICAL DATA SHEET

### InterMask™ : KN95 Face Mask, 12889-001, Protective Mask

<b>Total Mask</b>	This product consists of a mask body (made of non-woven fabric, hot air cotton and polypropylene melt-blown fabric) and a mask belt (elastic non-woven fabric), prevent dust, haze, direct penetration of pathogenic microorganisms, body fluids, and particulate matter provides a physical barrier.
-------------------	---

Packaging	Qty	Dimension	Weight
Primary Packing	2000 pcs / case		
Secondary Packing	50 Pcs in a box/polybag		

#### Instructions

1. Choose the appropriate mask size.
2. Perform hand hygiene before putting on a mask.
3. The mask should fit snugly over the face:
4. Most masks adopt a three-layer design. In general, the 2 coloured side/the side with folds facing downwards of the mask should face outwards with the metallic strip uppermost
5. For ear-loops type, position the elastic bands around both ears.
6. Extend the mask to fully cover mouth, nose and chin
7. Mould the metallic strip over nose bridge and the surgical mask should fit snugly over the face
8. Avoid touching the surgical mask after wearing. Otherwise, should perform hand hygiene before and after touching the mask.
9. For ear-loops type, hold both the ear loops and take-off gently from face. Avoid touching the outside of mask during taking-off as it may be covered with germs.
10. After taking off the surgical mask, discard in a lidded rubbish bin and perform hand hygiene immediately.
11. Change mask timely. In general, surgical mask should not be reused. Replace the mask immediately if it is damaged or soiled.

## TECHNICAL DATA SHEET

### InterMask™ : KN95 Face Mask, 12889-001, Protective Mask

**List of Item Codes and Attributes**

Item Code / Part Number	Description	Packaging Information	Color	Specific Attributes	Treatment
12889-001	InterMask™ : KN95 Protective Mask		White	Regular Earloop type	Sterile

**Disclaimer :**

This document and the information contained herein is the property of EcoPark (S) PTE. LTD. (Ecopark) and may not be used or reproduced without permission. Product users should conduct all appropriate testing or other evaluations to determine the suitability of Ecopark products for a particular purpose or use within a particular environment. ECOPARK DISCLAIMS ALL WARRANTIES OTHER THAN AS EXPRESSLY PROVIDED.