

# EcoSwab™ Disposable Virus Sampling Swab Kits

**P/N: 12899-002**

EcoSwab™ Disposable Virus Sampling Swab Kits includes sampling swabs that uses advanced jet implanted nylon fuzz technology and the sample can be easily penetrated, diluted and trapped in the fibre.

Country of Origin: China



## Technical Data Sheet

### General Information


Parameters	Specifications
Features	Nylon flocking Swab with outstanding sample collection and release capability. Release rate can reach 90%. Light Red Hank's Solution with antibiotics effectively inhibit reproduction of bacteria and fungi. Colorless and transparent preservative solutions effectively inactivate virus and reduce degradation of nucleic acid
Model	Nontoxic, Sterile
Swab Types	Nasal Swab, Pharyngeal Swab, Universal Swab

### Properties

Characteristics	Results
Tube Type	5ml/10ml
Specification	2ml/3ml
Rod Material	ABS
Tip Material	Flocked fabric
Dimension Overall Length	150±2mm (Pharyngeal, Nasal) 145±2mm(Universal)
Tip Diameter	Nasal - 3mm Pharyngeal - 5mm Universal - 3.5mm
Break off point	Nasal - 80mm Pharyngeal - 32mm Universal - 75 mm
Tip length	Nasal - 20mm Pharyngeal - 22mm Universal - 18 mm

### Packaging Details

Part #	Description	Packaging	NW/GW	Dimension
12899-002	EcoSwab™ Disposable Virus Sampling Swab Kits	Individually packaged	NA	NA



#### Disclaimer:

This document and the information contained herein is the property of EcoPark (S) PTE. LTD. (Ecopark) and may not be used or reproduced without permission. Product users should conduct all appropriate testing or other evaluations to determine the suitability of Ecopark products for a particular purpose or use within a particular environment. ECOPARK DISCLAIMS ALL WARRANTIES OTHER THAN WHAT HAS BEEN EXPRESSLY PROVIDED.

